

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Suppresses food intake, delays gastric emptying and decreases heat-induced edema. Might represent an endogenous ligand for maintaining homeostasis after stress.

Gene ID:

UCN2

Uniprot

Q96RP3

Synonyms:

Urocortin-2 (Stresscopin-related peptide) (Urocortin II) (Ucn II) (Urocortin-related peptide), UCN2, SRP URP


Immunogen:

Recombinant Human Urocortin-2 protein (72-112AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO38202 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO38202 at dilution of 1:100.