

# TTC5 Antibody, HRP conjugated



PACO38211

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Adapter protein involved in p53/TP53 response that acts by regulating and mediating the assembly of multi-protein complexes. Required to facilitate the interaction between JMY and p300/EP300 and increase p53/TP53-dependent transcription and apoptosis. Prevents p53/TP53 degradation by MDM2.

**Gene ID:**

TTC5

**Uniprot**

Q8N0Z6

**Synonyms:**

Tetratricopeptide repeat protein 5 (TPR repeat protein 5) (Stress-responsive activator of p300) (Strap), TTC5

**Immunogen:**

Recombinant Human Tetratricopeptide repeat protein 5 protein (201-440AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A