

TTC5 Antibody, FITC conjugated



PACO38212

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Adapter protein involved in p53/TP53 response that acts by regulating and mediating the assembly of multi-protein complexes. Required to facilitate the interaction between JMY and p300/EP300 and increase p53/TP53-dependent transcription and apoptosis. Prevents p53/TP53 degradation by MDM2.

Gene ID:

TTC5

Uniprot

Q8N0Z6

Synonyms:

Tetratricopeptide repeat protein 5 (TPR repeat protein 5) (Stress-responsive activator of p300) (Strap), TTC5

Immunogen:

Recombinant Human Tetratricopeptide repeat protein 5 protein (201-440AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A