

CPXCR1 Antibody, Biotin conjugated



PACO38221

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CPXCR1

Uniprot

Q8N123

Synonyms:

CPX chromosomal region candidate gene 1 protein (Cancer/testis antigen 77) (CT77), CPXCR1

Immunogen:

Recombinant Human CPX chromosomal region candidate gene 1 protein (1-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A