

C4orf33 Antibody




PACO38222

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse	Gene ID: C4orf33
Source: Rabbit	Uniprot Q8N1A6
Isotype: IgG	Synonyms: UPF0462 protein C4orf33, C4orf33
Applications: ELISA, WB	Immunogen: Recombinant Human UPF0462 protein C4orf33 protein (1-199AA).
Recommended dilutions: ELISA:1:2000-1:10000, WB:1:1000-1:5000	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: C4orf33 antibody at 4 μ g/ml. Lane 1: Mouse kidney tissue. Lane 2: Mouse liver tissue. Lane 3: Mouse lung tissue. Lane 4: Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 24 kDa. Observed band size: 24 kDa.