

# OSR2 Antibody, FITC conjugated



PACO38236

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

nucleus, sequence-specific DNA binding, bone morphogenesis, cell differentiation, chondrocyte differentiation, embryo development, embryonic digit morphogenesis, embryonic forelimb morphogenesis, embryonic hindlimb morphogenesis, embryonic skeletal joint development.

**Gene ID:**

OSR2

**Uniprot**

Q8N2R0

**Synonyms:**

Protein odd-skipped-related 2, OSR2

**Immunogen:**

Recombinant Human Protein odd-skipped-related 2 protein (101-276AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A