

OSR2 Antibody, Biotin conjugated



PACO38237

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

nucleus, sequence-specific DNA binding, bone morphogenesis, cell differentiation, chondrocyte differentiation, embryo development, embryonic digit morphogenesis, embryonic forelimb morphogenesis, embryonic hindlimb morphogenesis, embryonic skeletal joint development.

Gene ID:

OSR2

Uniprot

Q8N2R0

Synonyms:

Protein odd-skipped-related 2, OSR2

Immunogen:

Recombinant Human Protein odd-skipped-related 2 protein (101-276AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A