

PPM1K Antibody, Biotin conjugated



PACO38245

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulates the mitochondrial permeability transition pore and is essential for cellular survival and development.

Gene ID:

PPM1K

Uniprot

Q8N3J5

Synonyms:

Protein phosphatase 1K, mitochondrial (EC 3.1.3.16) (PP2C domain-containing protein phosphatase 1K) (PP2C-like mitochondrial protein) (PP2C-type mitochondrial phosphoprotein phosphatase) (PTMP) (Protein phosphatase 2C isoform kappa) (PP2C-kappa), PPM1K, PP2CM

Immunogen:

Recombinant Human Protein phosphatase 1K, mitochondrial protein (30-372AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A