

LIX1 Antibody, Biotin conjugated



PACO38253

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasm, autophagosome maturation, autophagy.

Gene ID:

LIX1

Uniprot

Q8N485

Synonyms:

Protein limb expression 1 homolog, LIX1, C5orf11

Immunogen:

Recombinant Human Protein limb expression 1 homolog protein (1-282AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A