

# WDTC1 Antibody, HRP conjugated



PACO38263

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.

**Gene ID:**

WDTC1

**Uniprot**

Q8N5D0

**Synonyms:**

WD and tetratricopeptide repeats protein 1, WDTC1, KIAA1037

**Immunogen:**

Recombinant Human WD and tetratricopeptide repeats protein 1 protein (201-420AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A