

WDTC1 Antibody, Biotin conjugated



PACO38265

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.

Gene ID:

WDTC1

Uniprot

Q8N5D0

Synonyms:

WD and tetratricopeptide repeats protein 1, WDTC1, KIAA1037

Immunogen:

Recombinant Human WD and tetratricopeptide repeats protein 1 protein (201-420AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A