

# PCDH20 Antibody, Biotin conjugated



PACO38281

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Potential calcium-dependent cell-adhesion protein.

**Gene ID:**

PCDH20

**Uniprot**

Q8N6Y1

**Synonyms:**

Protocadherin-20 (Protocadherin-13), PCDH20, PCDH13

**Immunogen:**

Recombinant Human Protocadherin-20 protein (61-360AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A