

PACO38282

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

EFHB

Uniprot

Q8N7U6

Synonyms:

EF-hand domain-containing family member B, EFHB


Immunogen:

Recombinant Human EF-hand domain-containing family member B protein (405-644AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: EFHB antibody at 4 μ g/ml. Lane 1: HepG2 whole cell lysate. Lane 2: PC-3 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 94, 71, 80 kDa. Observed band size: 94 kDa.