

TRIML1 Antibody, FITC conjugated



PACO38296

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable E3 ubiquitin-protein ligase which plays an important role in blastocyst development.

Gene ID:

TRIML1

Uniprot

Q8N9V2

Synonyms:

Probable E3 ubiquitin-protein ligase TRIML1 (EC 2.3.2.27) (RING finger protein 209) (RING-type E3 ubiquitin transferase TRIML1) (Tripartite motif family-like protein 1), TRIML1, RNF209

Immunogen:

Recombinant Human Probable E3 ubiquitin-protein ligase TRIML1 protein (277-468AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A