

WFDC9 Antibody



PACO38326

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

WFDC9

Uniprot

Q8NEX5

Synonyms:

Protein WFDC9, WFDC9, WAP9

Immunogen:

Recombinant Human Protein WFDC9 protein (24-89AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A