

# GSDMB Antibody, FITC conjugated



PACO38340

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role as secretory or metabolic product involved in secretory pathway. May also play a role in achieving and maintaining the final differentiation state of epithelial cells.

**Gene ID:**

GSDMB

**Uniprot**

Q8TAX9

**Synonyms:**

Gasdermin-B (Gasdermin-like protein), GSDMB, GSDML

**Immunogen:**

Recombinant Human Gasdermin-B protein (63-310AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A