

GSDMB Antibody, Biotin conjugated



PACO38341

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role as secretory or metabolic product involved in secretory pathway. May also play a role in achieving and maintaining the final differentiation state of epithelial cells.

Gene ID:

GSDMB

Uniprot

Q8TAX9

Synonyms:

Gasdermin-B (Gasdermin-like protein), GSDMB, GSDML

Immunogen:

Recombinant Human Gasdermin-B protein (63-310AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A