

# DENND1A Antibody



PACO38354

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Guanine nucleotide exchange factor (GEF) regulating clathrin-mediated endocytosis through RAB35 activation. Promotes the exchange of GDP to GTP, converting inactive GDP-bound RAB35 into its active GTP-bound form. Regulates clathrin-mediated endocytosis of synaptic vesicles.

**Gene ID:**

DENND1A

**Uniprot**

Q8TEH3

**Synonyms:**

DENN domain-containing protein 1A (Connecdenn 1) (Connecdenn) (Protein FAM31A), DENND1A, FAM31A KIAA1608

**Immunogen:**


Recombinant Human DENN domain-containing protein 1A protein (1-559AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO38354 at dilution of 1:100.