

PTCD2 Antibody, HRP conjugated



PACO38371

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in mitochondrial RNA maturation and mitochondrial respiratory chain function.

Gene ID:

PTCD2

Uniprot

Q8WV60

Synonyms:

Pentatricopeptide repeat-containing protein 2, mitochondrial, PTCD2

Immunogen:

Recombinant Human Pentatricopeptide repeat-containing protein 2, mitochondrial protein (153-388AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A