

# PTCD2 Antibody, Biotin conjugated



PACO38373

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in mitochondrial RNA maturation and mitochondrial respiratory chain function.

**Gene ID:**

PTCD2

**Uniprot**

Q8WV60

**Synonyms:**

Pentatricopeptide repeat-containing protein 2, mitochondrial, PTCD2

**Immunogen:**

Recombinant Human Pentatricopeptide repeat-containing protein 2, mitochondrial protein (153-388AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A