

Serpina3n Antibody, HRP conjugated



PACO38387

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

blood microparticle, extracellular exosome, extracellular region, extracellular space, serine-type endopeptidase inhibitor activity, negative regulation of endopeptidase activity, response to cytokine, response to peptide hormone.

Gene ID:

Serpina3n

Uniprot

Q91WP6

Synonyms:

Serine protease inhibitor A3N (Serp A3N), Serpina3n, Spi2

Immunogen:

Recombinant Mouse Serine protease inhibitor A3N protein (21-418AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A