

ZKSCAN4 Antibody



PACO38398

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the transcriptional activation of MDM2 and EP300 genes.

Gene ID:

ZKSCAN4

Uniprot

Q969J2

Synonyms:

Zinc finger protein with KRAB and SCAN domains 4 (P373c6.1) (Zinc finger protein 307) (Zinc finger protein 427), ZKSCAN4, ZNF307 ZNF427

Immunogen:

Recombinant Human Zinc finger protein with KRAB and SCAN domains 4 protein (264-545AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A