

CLUAP1 Antibody, HRP conjugated



PACO38415

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for cilia biogenesis. Appears to function within the multiple intraflagellar transport complex B (IFT-B). Key regulator of hedgehog signaling.

Gene ID:

CLUAP1

Uniprot

Q96AJ1

Synonyms:

Clusterin-associated protein 1 (Qilin), CLUAP1, KIAA0643

Immunogen:

Recombinant Human Clusterin-associated protein 1 protein (161-413AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A