

ZMYND19 Antibody



PACO38434

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved as a regulatory molecule in GPR24/MCH-R1 signaling.

Gene ID:

ZMYND19

Uniprot

Q96E35

Synonyms:

Zinc finger MYND domain-containing protein 19 (Melanin-concentrating hormone receptor 1-interacting zinc finger protein) (MCH-R1-interacting zinc finger protein), ZMYND19, MIZIP


Immunogen:

Recombinant Human Zinc finger MYND domain-containing protein 19 protein (1-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO38434 at dilution of 1:100.