

DDX27 Antibody



PACO38462

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable ATP-dependent RNA helicase.

Gene ID:

DDX27

Uniprot

Q96GQ7

Synonyms:

Probable ATP-dependent RNA helicase DDX27 (EC 3.6.4.13) (DEAD box protein 27), DDX27, cPERP-F RHLF

Immunogen:

Recombinant Human Probable ATP-dependent RNA helicase DDX27 protein (30-794AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A