

# DDX27 Antibody, FITC conjugated



PACO38464

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probable ATP-dependent RNA helicase.

**Gene ID:**

DDX27

**Uniprot**

Q96GQ7

**Synonyms:**

Probable ATP-dependent RNA helicase DDX27 (EC 3.6.4.13) (DEAD box protein 27), DDX27, cPERP-F RHLF

**Immunogen:**

Recombinant Human Probable ATP-dependent RNA helicase DDX27 protein (30-794AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A