

PRR11 Antibody, HRP conjugated



PACO38475

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a critical role in cell cycle progression.

Gene ID:

PRR11

Uniprot

Q96HE9

Synonyms:

Proline-rich protein 11, PRR11

Immunogen:

Recombinant Human Proline-rich protein 11 protein (1-360AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A