

# PRR11 Antibody, Biotin conjugated



PACO38477

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays a critical role in cell cycle progression.

**Gene ID:**

PRR11

**Uniprot**

Q96HE9

**Synonyms:**

Proline-rich protein 11, PRR11

**Immunogen:**

Recombinant Human Proline-rich protein 11 protein (1-360AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A