

FAM167A Antibody



PACO38482

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

FAM167A

Uniprot

Q96KS9

Synonyms:

Protein FAM167A, FAM167A, C8orf13


Immunogen:

Recombinant Human Protein FAM167A protein (1-214AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO38482 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO38482 at dilution of 1:100.