

# ART5 Antibody, Biotin conjugated



PACO38489

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ART5

**Uniprot**

Q96L15

**Synonyms:**

Ecto-ADP-ribosyltransferase 5 (EC 2.4.2.31) (ADP-ribosyltransferase C2 and C3 toxin-like 5) (ARTC5) (Mono(ADP-ribosyl)transferase 5) (NAD(P)(+)-arginine ADP-ribosyltransferase 5), ART5

**Immunogen:**

Recombinant Human Ecto-ADP-ribosyltransferase 5 protein (45-291AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A