

# UBXN10 Antibody



PACO38490

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

UBXN10

**Uniprot**

Q96LJ8

**Synonyms:**

UBX domain-containing protein 10 (UBX domain-containing protein 3), UBXN10, UBXD3

**Immunogen:**


Recombinant Human UBX domain-containing protein 10 protein (1-280AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO38490 at dilution of 1:100.