

FAM19A4 Antibody



PACO38494

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FAM19A4

Uniprot

Q96LR4

Synonyms:

Protein FAM19A4 (Chemokine-like protein Tafa-4), FAM19A4, Tafa4

Immunogen:

Recombinant Human Protein FAM19A4 protein (36-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A