

FCRL4 Antibody, FITC conjugated



PACO38504

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May function as an inhibitor of the B-cell receptor signaling. May function in the B-cell-mediated immune response.

Gene ID:

FCRL4

Uniprot

Q96PJ5

Synonyms:

Fc receptor-like protein 4 (FcR-like protein 4) (FcRL4) (Fc receptor homolog 4) (FcRH4) (IFGP family protein 2) (hIFGP2) (Immune receptor translocation-associated protein 1) (CD antigen CD307d), FCRL4, FCRH4 IFGP2 IRTA1

Immunogen:

Recombinant Human Fc receptor-like protein 4 protein (20-387AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A