

# FCRL4 Antibody, Biotin conjugated



PACO38505

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May function as an inhibitor of the B-cell receptor signaling. May function in the B-cell-mediated immune response.

**Gene ID:**

FCRL4

**Uniprot**

Q96PJ5

**Synonyms:**

Fc receptor-like protein 4 (FcR-like protein 4) (FcRL4) (Fc receptor homolog 4) (FcRH4) (IFGP family protein 2) (hIFGP2) (Immune receptor translocation-associated protein 1) (CD antigen CD307d), FCRL4, FCRH4 IFGP2 IRTA1

**Immunogen:**

Recombinant Human Fc receptor-like protein 4 protein (20-387AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A