

LYSMD1 Antibody



PACO38506

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

LYSMD1

Uniprot

Q96S90

Synonyms:

LysM and putative peptidoglycan-binding domain-containing protein 1, LYSMD1


Immunogen:

Recombinant Human LysM and putative peptidoglycan-binding domain-containing protein 1 protein (1-227AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO38506 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO38506 at dilution of 1:100.