

# C8orf59 Antibody



PACO38514

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

mitophagy in response to mitochondrial depolarization, positive regulation of defense response to virus by host, xenophagy.

**Gene ID:**

C8orf59

**Uniprot**

Q8N0T1

**Synonyms:**

Uncharacterized protein C8orf59, C8orf59

**Immunogen:**

Recombinant Human Uncharacterized protein C8orf59 protein (1-43AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A