

PACO38518

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Enhances the solubility and secretion of VASH1 and prevents its ubiquitination. Also enhances secretion of VASH2.

Gene ID:

SVBP

Uniprot

Q8N300

Synonyms:

Small vasohibin-binding protein (Coiled coil domain-containing protein 23), SVBP, CCDC23

Immunogen:

Recombinant Human Small vasohibin-binding protein (1-66AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A