

# HDDC3 Antibody, HRP conjugated



PACO38539

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

ppGpp hydrolyzing enzyme involved in starvation response.

**Gene ID:**

HDDC3

**Uniprot**

Q8N4P3

**Synonyms:**

Guanosine-3',5'-bis(diphosphate) 3'-pyrophosphohydrolase MESH1 (EC 3.1.7.2) (HD domain-containing protein 3) (Metazoan SpoT homolog 1) (MESH1) (Penta-phosphate guanosine-3'-pyrophosphohydrolase) ((ppGpp)ase), HDDC3, MESH1

**Immunogen:**

Recombinant Human Guanosine-3',5'-bis(diphosphate) 3'-pyrophosphohydrolase MESH1 protein (2-140AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A