

HDDC3 Antibody, Biotin conjugated



PACO38541

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

ppGpp hydrolyzing enzyme involved in starvation response.

Gene ID:

HDDC3

Uniprot

Q8N4P3

Synonyms:

Guanosine-3',5'-bis(diphosphate) 3'-pyrophosphohydrolase MESH1 (EC 3.1.7.2) (HD domain-containing protein 3) (Metazoan SpoT homolog 1) (MESH1) (Penta-phosphate guanosine-3'-pyrophosphohydrolase) ((ppGpp)ase), HDDC3, MESH1

Immunogen:

Recombinant Human Guanosine-3',5'-bis(diphosphate) 3'-pyrophosphohydrolase MESH1 protein (2-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A