

VGLL2 Antibody



PACO38558

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May act as a specific coactivator for the mammalian TEFs. May play a role in the development of skeletal muscles.

Gene ID:

VGLL2

Uniprot

Q8N8G2

Synonyms:

Transcription cofactor vestigial-like protein 2 (Vgl-2) (Protein VITO1), VGLL2, VITO1


Immunogen:

Recombinant Human Transcription cofactor vestigial-like protein 2 protein (1-143AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO38558 at dilution of 1:100.