

L3MBTL4 Antibody, FITC conjugated



PACO38564

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Putative Polycomb group (PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility.

Gene ID:

L3MBTL4

Uniprot

Q8NA19

Synonyms:

Lethal(3)malignant brain tumor-like protein 4 (H-l(3)mbt-like protein 4) (L(3)mbt-like protein 4) (L3mbt-like 4), L3MBTL4

Immunogen:

Recombinant Human Lethal(3)malignant brain tumor-like protein 4 protein (1-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A