

# EHBP1 Antibody



PACO38570

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in actin reorganization. Links clathrin-mediated endocytosis to the actin cytoskeleton.

**Gene ID:**

EHBP1

**Uniprot**

Q8NDI1

**Synonyms:**

EH domain-binding protein 1, EHBP1, KIAA0903 NACSIN

**Immunogen:**

Recombinant Human EH domain-binding protein 1 protein (1-106AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A