

EHBP1 Antibody, HRP conjugated



PACO38571

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in actin reorganization. Links clathrin-mediated endocytosis to the actin cytoskeleton.

Gene ID:

EHBP1

Uniprot

Q8NDI1

Synonyms:

EH domain-binding protein 1, EHBP1, KIAA0903 NACSIN

Immunogen:

Recombinant Human EH domain-binding protein 1 protein (1-106AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A