

PACO38574

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the biosynthesis of ganglioside GD1A from GM1B. Transfers CMP-NeuAc with an alpha-2,6-linkage to GalNAc residue on NeuAc-alpha-2,3-Gal-beta-1,3-GalNAc of glycoproteins and glycolipids. ST6GalNAcIII prefers glycolipids to glycoproteins.

Gene ID:

ST6GALNAC3

Uniprot

Q8NDV1

Synonyms:

Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 3 (EC 2.4.99. -) (GalNAc alpha-2,6-sialyltransferase III) (ST6GalNAc III) (ST6GalNAcIII) (STY) (Sialyltransferase 7C) (SIAT7-C), ST6GALNAC3, SIAT7C

Immunogen:

Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 3 protein (40-305AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A