

PPTC7 Antibody



PACO38594

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

PPTC7

Uniprot

Q8NI37

Synonyms:

Protein phosphatase PTC7 homolog (EC 3.1.3.16) (T-cell activation protein phosphatase 2C) (TA-PP2C) (T-cell activation protein phosphatase 2C-like), PPTC7, TAPP2C


Immunogen:

Recombinant Human Protein phosphatase PTC7 homolog protein (1-304AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO38594 at dilution of 1:100.