

SCGB1C1 Antibody



PACO38602

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SCGB1C1

Uniprot

Q8TD33

Synonyms:

Secretoglobin family 1C member 1 (Secretoglobin RYD5), SCGB1C1

Immunogen:

Recombinant Human Secretoglobin family 1C member 1 protein (24-95AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A