

CXorf40A Antibody



PACO38614

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May have an important role of cell protection in inflammation reaction.

Gene ID:

CXorf40A

Uniprot

Q8TE69

Synonyms:

Protein CXorf40A (Endothelial-overexpressed lipopolysaccharide-associated factor 1), CXorf40A, CXorf40 EOLA1

Immunogen:

Recombinant Human Protein CXorf40A protein (21-158AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A