

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

The RIC1-RGP1 complex acts as a guanine nucleotide exchange factor (GEF), which activates RAB6A by exchanging bound GDP for free GTP and may thereby be required for efficient fusion of endosome-derived vesicles with the Golgi compartment. The RIC1-RGP1 complex participates in the recycling of mannose-6-phosphate receptors.

**Gene ID:**

RGP1

**Uniprot**

Q92546

**Synonyms:**

RAB6A-GEF complex partner protein 2 (Retrograde Golgi transport protein RGP1 homolog), RGP1, KIAA0258

**Immunogen:**

Recombinant Human RAB6A-GEF complex partner protein 2 protein (56-350AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A