

# RPGR Antibody, FITC conjugated



PACO38668

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Could be a guanine-nucleotide releasing factor. Plays a role in ciliogenesis. Probably regulates cilia formation by regulating actin stress filaments and cell contractility. Plays an important role in photoreceptor integrity. May play a critical role in spermatogenesis and in intraflagellar transport processes. May be involved in microtubule organization and regulation of transport in primary cilia.

**Gene ID:**

RPGR

**Uniprot**

Q92834

**Synonyms:**

X-linked retinitis pigmentosa GTPase regulator, RPGR, RP3 XLRP3

**Immunogen:**

Recombinant Human X-linked retinitis pigmentosa GTPase regulator protein (1-180AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A