

ARRDC3 Antibody, Biotin conjugated



PACO38693

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Adapter protein that plays a role in regulating cell-surface expression of adrenergic receptors and probably also other G protein-coupled receptors. Plays a role in NEDD4-mediated ubiquitination and endocytosis of activated ADRB2 and subsequent ADRB2 degradation. May recruit NEDD4 to ADRB2. Alternatively, may function as adapter protein that does not play a major role in recruiting NEDD4 to ADRB2, but rather plays a role in targeting ADRB2 to endosomes.

Gene ID:

ARRDC3

Uniprot

Q96B67

Synonyms:

Arrestin domain-containing protein 3 (TBP-2-like inducible membrane protein) (TLIMP), ARRDC3, KIAA1376

Immunogen:

Recombinant Human Arrestin domain-containing protein 3 protein (1-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A